



Enantiomeric Separation of Isoprenaline with HP-β-Cyclodextrins as Chiral Selector



Capillary Electrophoresis is a powerful analytical technique for chiral analysis.



Application developed by the Molecular Pharmacology Department, Grenoble-Alpes University
corinne.ravelet@univ-grenoble-alpes.fr


INTRODUCTION

Racemic Isoprenaline is a catecholamine used for the treatment of bradycardia and heart block. Chiral separation is mandatory for the quantitation of each enantiomer. In this application note, we use the Wyn-CE Capillary Electrophoresis system, for the isoprenaline chiral separation. Effects of pH buffer, temperature and cyclodextrin concentrations on the separation are described.

STANDARD AND REAL ANALYSIS

Buffer : Phosphate buffer + HP-β-CD, variable pH
Capillary : bubble bare-fused silica, L = 64,5 cm, l = 56 cm, ID = 50 μm
Injection : hydrodynamic, 50 mbar, 4 s
Voltage : +30 kV
Detection : UV, 200 nm
Temperature : 35 °C
Sample concentration : 250 mg/L


Principle



Needs for Chiral Separation :


- 1) Different complex stabilities for each enantiomer
- 2) Different mobilities between free and complexed species

The used CD



Hydroxypropyl-β-cyclodextrin (HP-β-CD)

Effect of HP-β-CD concentration



Effect of pH

